

# SAS™ RSVAlert Test

CLIA  
WAIVED



► **CPT Code:**  
87807/87807QW\*

► **Sensitivity:**  
97%

► **Specificity:**  
98%

► **Specimen:**  
Nasal Pharyngeal  
Washes/Aspirates  
& Swabs

SAS™ RSVAlert Test is a rapid lateral-flow, membrane based immunoassay for the qualitative detection of Respiratory Syncytial Virus. RSV is the most common cause of lower respiratory tract infections such as bronchiolitis and pneumonia among infants and young children. RSV can cause serious lower respiratory tract infections throughout life, especially among those with compromised respiratory, cardiac or immune systems.



- CLIA Classification: Waived
- Simple, 15-minute rapid, visual, membrane test
- Optional extraction step for bloody or mucoidal specimens
- Room temperature storage
- 18-month shelf life

\*"QW" modifier is added to report use of CLIA-waived test systems for Medicare/Medicaid claims.